





The latest evolution of the unmistakable and historical Linda 3 model has arrived: LINDA EVO is finally here. **LINDA EVO** is the ultimate mid-range knee-break chair, perfectly suitable to treat all patients, even children, elderly people, and persons with reduced mobility, for both dental and medical procedures. The distinctive features of the well renown Tecnodent kneebreak chair concept, such as compactness, easy accessibility, clean lines, strong solidity, and comfort of use, are preserved and empowered with modern design, technology, and structure. Low voltage motors, with Hall sensor limit switches, guarantee smoothness and perfect movement fluidity, thanks also to softstart and soft-stop features, while the unmistakable designer armrests and split-section upholstery, along with the enveloping backrest, uphold the smart design statement.

Versatile and comfortable, **LINDA EVO** embodies the perfect "multifunctional smart knee-break", thanks to its elegant yet

compact, ergonomic, and functional style.

LINDA EVO STANDARD FEATURES

- Fully electromechanical upward, downward and tilt movements of the backrest
- Backrest movement compensated, up to horizontal position
- Tilting of the backrest synchronized with the elongation of the leg rest
- Double-articulated headrest model VALE
- High enveloping backrest
- Complete symmetry, usable by both right and left-handed professionals
- Safety stops on backrest and lifting arm

- Automatic program for complete return to the patient lowering position ("zero")
- Three user programmable customized working positions
- Automatic rinse program with return to the last working position
- Low voltage controls: 5 V
- Joystick control, on the base
- Soft-start and soft-stop movement



LINDA EVO - Options

- Double-articulated headrest model ELLE MOBILE
- Head rest model ELLE-2 MOBILE
- Magnetic cushion C95
- Magnetic cushion C2002
- Multi-articulated headrest model "UNI tappezzato"
- Multi-articulated headrest model "LIGHT"
- · Right and/or left swiveling armrest
- Electronic board with lamp activation relay
- Rotation system around the vertical axis under the chair base (320°)
- Infrared hand remote control
- Children booster seat
- ENT footrest
- Unit / lamp bracket
- Seamless memory foam upholstery
- Handmade soft-foam upholstery
- Handmade soft-foam upholstery with contrasting color stitching



LINDA EVO SUR and AES set-ups

- Lateral tool holder bars
- Adjustable surgical rotating armrests
- COMFORT wide multi-articulated headrest
- PRIMA FLEX compact directional lamp
- Polyurethane gel handrest
- Paper roll holder
- Tablet / Articulating Mirror holder



Joystick pedal control



Soft Foam

Swiveling armrests





Base with Rotation system

IR remote control





ENT footrest

STANDARD AND OPTIONAL CHARACTERISTICS		NEXT	A EVO
• Standard • Optional		LINDA NEXT	LINDA EVC
Headrests	Vale	•	•
	ELLE/ELLE2 MOBILE	•	•
	UNI/Light	•	•
	COMFORT	•	•
Pillows	C95	•	•
	C2002	•	•
Armrest Backrest Pillows	Backrest motorized height adjustment	•	
	Large enveloping	•	
	High enveloping		•
Armrest	Right armrest	•	•
	Left armrest	•	•
Programs	Automatic built in programs: Patient	•	
	Positionning, Table, Trendelenburg		
	9 user programmable positions (3x profile)	•	
	3 user programmable positions		•
	Return to "zero" program	•	•
	Rinse program	•	•
Other features	Rotation	•	•
	Lamp / unit bracket	•	•
	Back-up Battery	•	
	Wheels	•	
Function	Trendelenburg	•	•
	Knee-break legrest	•	•
	Independent legrest movement	•	
	Securities	•	•
Control	Joystick	•	•
	Infrared remote control	•	•

HEADRESTS

Tecnodent manual Headrests are available in 6 different models: ELLE/ELLE2 MOBILE, UNI, Light, COMFORT and VALE. Each headrest was made for meeting every operator's need.

ELLE/ELLE2 MOBILE is an anatomical headrest, characterized by movement of inclination of the upper part (double articulation) for a correct cervical support.



ELLE/ELLE2





C2002

Magnetic pillows. ELLE/ELLE2 headrests can be combined with magnetic pillow C95 or C2002. Either pillow can be applied directly on the backrest, for ensure a better comfort, independently of the patient's height.







COMFORT

Uni, COMFORT and Light. The UNI, COMFORT and Light headrests are multiarticulated and allow positioning of the headrest along 3 axis.



VALE is a double articulated headrest with enveloping cushion.

VALE

