





There is a brand new star in Tecnodent family... Eco Lux.

Eco Lux is a perfect match between contemporary shapes, solid structure and charming Italian design.

The gentle line of the upper part of the chair with its design armrests and elegant structural shining polyurethane meets the renewed unmistakable long lasting lower pantograph arm and base made with integral PU.

Eco Lux embodies the essence of The Dental Chair combining a unique ergonomic with a perfect comfortable experience for both doctor and patient.

Eco Lux

- Anatomical headrest ELLE2 MOBILE in integral PU with mouting along vertical axis and double articulation for a patient's comfortable posture and a perfect view during the treatment;
- Automatic programs: "zero", rinse, last position;
- Compensated backrest tilting, up to the horizontal position;
- Complete symmetry for right-handed or left-handed professionals;
- Electromechanical lifting, lowering and tilting movements;
- Joystick control on the base for fast and efficient movements of the chair without using the hands, assuring a perfect hygiene;
- Large Integrated utility box;
- Narrow and slim anatomical backrest for the maximum freedom of action of the professional;
- Possible customization with different colors of polyurethane and leatherette fabric;
- Security STOPS on the backrest, on the base and on the pantograph arm;
- Wrapping seamless upholstery in high quality leatherette fabric.

OPTIONS

- Children pillow available;
- Different headrests are available according to specific treatment needs;
- European or American unit bracket;
- IR remote control;
- Large wrapping backrest;
- Memory Foam and Hand Made Soft Foam upholstery;
- Mobile foot control;
- Possibility to program 3 different working positions;
- Possibility to equip the chair with one/two rotating armrests to improve the accessibility of the chair and to assure a natural support for the arms, allowing a relaxed position during the treatment;
- Rotation and Front adaptor;
- Right drilling/ rear drilling/ double side drilling/ no drilling;
- Trendelenburg position.



Narrow backrest

This backrest is studied to allow the doctor the closest approach to the patient during the treatment, also guaranteeing an excellent comfort to the patient.

Mobile foot control "Easy Touch" (OPTION)

Easy and functional foot control which provides different precise treatment positions and doctor's hands-free control.

Infrared remote control (OPTION)

Manual control of the movements and set to reach programmed positions.



Large anatomical backrest (OPTION)

Guarantees the maximum comfort, wrapping the patient's back in a pleasant relax.

Pedal control

A high precision joystick pedal controls directly in a simple and efficient way the movements without using the hands.

Right and/or left Lux it armrest (OPTION)

Rotating armrests help the patient to get on/off the chair.

TECHNICAL FEATURES

| Maximum height | mm | Without rotation 780 | With rotation 795 | | | | |
|---|----|-------------------------|----------------------|--|--|--|--|
| Minimum height | mm | Without rotation 380 | With rotation 395 | | | | |
| Backrest tilting angle | 0 | -15/+65 | | | | | |
| Rotation angle | 0 | ±30 | | | | | |
| Width | mm | 660 | | | | | |
| Length (without headrest) | mm | 1570 | | | | | |
| Length (with headrest) | mm | 2050 | | | | | |
| Gross weight | Kg | 170 | | | | | |
| Net weight | Kg | 140 | | | | | |
| Maximum lifting capability (patient + unit) | Kg | 135+75 | | | | | |
| Classification | | 1 | | | | | |
| Туре | | В | | | | | |



STANDARD AND OPTIONAL CHARACTERISTICS

| | • Standard • Optional | 2009next | Eco Lux | ECO.next | Linda nex | Sting |
|-------------------|--|----------|---------|----------|-----------|-------|
| Pillows Headrests | ELLE FISSO | • | | | • | |
| | ELLE/ELLE2 MOBILE | • | • | • | • | • |
| | UNI/Light | • | • | • | • | • |
| Pillows | C95 | • | • | • | • | • |
| | C2002 | • | • | • | • | • |
| Armrest Backrest | Backrest motorized height adjustment | | | | • | |
| | Narrow | | • | • | | |
| | Large | • | • | • | • | • |
| | Right armrest | • | • | • | • | • |
| | Left armrest | • | • | • | • | • |
| | 3 customizable programs (T2000 PCB) | • | • | • | | • |
| Programs | Automatic built in programs: Patient | | | | • | |
| | Positionning, Table, Trendelenburg | | | | | |
| | 9 user programmable positions (3x profile) | | | | • | |
| | Return to "zero" program | • | • | • | • | • |
| | Rinse program | • | • | • | • | • |
| | Double side drilled upper body | • | • | • | | • |
| tion | Fixed central arm | | • | • | | |
| Unit connection | Rotation | • | • | • | • | • |
| | Rotation with central arm | • | • | • | | • |
| | Rear connection for unit bracket | • | • | • | • | • |
| | Trendelenburg | • | • | • | • | • |
| Function | Knee-break legrest | | | | • | |
| | Independent legrest movement | | | | • | |
| | Securities | • | • | • | • | • |
| | Small integrated utility box | | | | | |
| | Large integrated utility box | • | • | • | | • |
| Control | Joystick | • | • | • | • | • |
| | Mobile foot control | • | • | • | • | • |
| | Infrared remote control | | • | • | • | • |

HEADRESTS

Tecnodent manual Headrests are available in 4 different models: ELLE FISSO, ELLE/ELLE2 MOBILE, UNI and Light. Every headrest was made for meeting every operator's need.

ELLE/ELLE2. ELLE/ELLE2 FISSO has only the movment along its vertical axis. ELLE/ELLE2 MO-BILE is an anatomical headrest, characterized by movement of inclination of the upper part (double articulation) for a correct cervical support.



ELLE/ELLE2





Magnetic pillows. ELLE/ELLE2 type, both fix and mobile version, can be combined with magnetic pillow C95 or C2002. Either pillow can be applied directly on the backrest, for ensure a better comfort, independently of the patient's height.





UNI

UNI and Light. The UNI and Light headrests are multi-articulated and allows positioning of the headrest along 3 axis.

