







EGG R stool

- High-resistant polyurethane base equipped with 5 rubber casters (also available with friction);
- Ergonomic and seamless round seat upholstery, with high quality leatherette fabric;
- Optional polyurethane footrest ring;
- The revolving and folding "U"shaped backrest easily turns into a comfortable armrest;
- Height adjustable with gas pump;
- Casters spacers can be mounted on the base to increase stool height;
- Possibility to customize the stool in different polyurethane colors.



EGG Iron stool

- Chromed aluminium design base equipped with5 rubber casters (also available with friction);
- Ergonomic and seamless round seat upholstery, with high quality leatherette fabric;
- Optional chromed aluminium footrest ring;
- The revolving and folding "U"-shaped backrest easily turns into a comfortable armrest;
- Height adjustable with gas pump, two different extensions available: standard (minimum) 564 mm and extra (maximum) 696;
- Possibility to customize the stool in different polyurethane colors;
- EGG can be also available in a "surgical" version, equipped with antistatic casters, fireproof and seamless upholsteries, up&down movement with hydraulic pedal pump and locking device to fasten the stool to the ground.



Saddle StingLE

- Chromed aluminium design base equipped with 5 rubber casters (also available with friction);
- Anatomical saddle seat with adjustable height and tilting control;
- Height adjustable with gas pump, two different extensions available: standard (minimum) 564 mm and extra (maximum) 696;
- Seamless upholstery with high quality leatherette fabric;
- Optional chromed aluminium footrest ring;
- Tilting rigidity control according to the body weight;
- Sting stools can be also available in a "surgical" version, equipped with antistatic casters, fireproof and seamless upholsteries, up&down movement with hydraulic pedal pump and locking device to fasten the stool to the ground.



StingST stool

- Chromed aluminium design base equipped with 5 rubber casters (also available with friction);
- Anatomical saddle seat with adjustable height and tilting control;
- Height adjustable with gas pump, two different extensions available: standard (minimum) 564 mm and extra (maximum) 696;
- Height-adjustable and anatomical backrest for perfect dorsal support;
- Seamless upholstery with high quality leatherette fabric;
- Optional chromed aluminium footrest ring;
- Tilting rigidity control according to the body weight;
- Sting stools can be also available in a "surgical" version, equipped with antistatic casters, fireproof and seamless upholsteries, up&down movement with hydraulic pedal pump and locking device to fasten the stool to the ground.



CHARACTERISTICS		GEL Line				EGG Line		ERGO Line		
• Standard • Optional	Pluto	PlutoLE	Venus	StingST	EGG R	EGG Iron	D	А	S	
Castors without brake	•	•	•	•	•	•	•	•	•	
Castors with brake	•	•	•	•	•	•	•	•	•	
Special polyurethane base	•	•	•							
Polyurethane "umbrella" base					•					
Aluminium base			•	•		•	•	•	•	
Footrest ring	•	•	•	•	•	•	•	•	•	
Spacers					•					
Gas spring of variable heights			•	•		•	•	•	•	
Round seat			•		•	•		•		
Saddle seat	•	•		•					•	
Anatomical backrest	•						•			
"U" shaped backrest			•		•	•		•		
Backrest adjustable in height	•			•			•			
Backrest inclination independent of seat										
inclination										
Inclination weight control	•	•		•						

TECHNICAL FEATURES		EGG R	EGG Iron	StingST
Minimum height	mm	470	510	530
Maximum height	mm	650	650	720
External castors diameter	mm	540	640	640
Gross weight	kg	25	18	22
Net weight	kg	22	14	19
Seat width	mm	380	380	400x370
Backrest inclination adjustment	0	-	-	24
Seat inclination adjustment	0	-	-	-8 / +6
Maximum lifting weight	kg	135	135	135

